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RU

MediVital®

*Illustrative picture
MediVital® with guard*

*Illustrative picture
MediVital® without guard*

1. FOREWORD

This medical device classified as class IIb complies the General Safety and Performance Requirements of the Medical Device Regulation (EU) 2017/745 and also comply with Transportable Pressure Equipment Directive 2010/35/EU and CDG TPE (Amendment) & (EU Exit) Regulations 2020.

Pressure regulators integrated with cylinder valves (also referred to as “combination valves”) are designed according to EN ISO 10524-3 and EN ISO 10297 standards.

Note: The precautions/warnings specified for oxygen are valid for all gases (gas specification is stated in Appendix Nr. 1), if not specified differently.

2. INTENDED USE

Medical combination valves MediVital® are designed to be fitted to gas cylinders used for medical gases up to 300 bar. These combination valves together with gas cylinders form gas packages used as gas supply points for medical devices and allow management of the medical gases and their mixtures supply. Medical combination valves are used to reduce high cylinder pressure to lower pressure or to specific outlet flow suitable for the administration of medical gases in the treatment of patients.

MediVital® combination valve is not intended for use with flammable anaesthetics and substances.

3. OPERATIONAL, TRANSPORT AND STORAGE SAFETY REQUIREMENTS

Always maintain oxygen cleanliness standards:

KEEP THE MEDICAL DEVICE AND ITS ASSOCIATED EQUIPMENT AWAY FROM:

- Heat sources (fire, cigarettes...),
- Flammable materials,
- Oil or grease (take a great care in the use of hand creams),
- Water,
- Dust.

The medical device and its associated equipment must be prevented from tipping over, turning over or falling.

Use the medical device and its associated equipment in well ventilated area only. Otherwise, there is a danger of ignition, which can damage and dysfunction the medical device and endanger the user.

Before the first use the medical device shall be kept in its original package. If removed from service (for transport, storage) GCE recommends using the original package (including inner packing materials). National laws, rules and regulations for medical gases, accident prevention and environmental protection must be observed.

OPERATING CONDITIONS		STORAGE AND TRANSPORT CONDITIONS	
	-20* / +65 °C		-40 / +65 °C

*Note: For inner tightness of the shut-off valve, during transport and storage of the combination valve mounted on a cylinder, the valid lower temperature limit is -40°C.

 In case of the combination valve storage temperature drops below -20°C do not use the combination valve until its temperature reaches at least -20°C .

 For the combination valves designed to be used with mixture of gases $\text{O}_2+\text{N}_2\text{O}$, the lowest operating temperature is $+5^{\circ}\text{C}$. In normal use of the combination valve, frosting can appear on the combination valve surface. This is caused by the gas inside the combination valve cooling when high pressure gas is being reduced to a low pressure (Joule-Thomson effect). Check that all patient associated equipment connected to the combination valve is via a hose of at least 2 metres length.

 $\text{O}_2+\text{N}_2\text{O}$ mixtures are temperature sensitive. N_2O begins to separate out from the mixture if the temperature falls below approximately -6°C . A homogenous mixture is achieved if the temperature has raised above 10°C and the cylinder is agitated. Before use, to ensure it is properly mixed, cylinders should be stored horizontally for 24 hours at a temperature above 10°C . If this is not practicable, before use the cylinders must be maintained at a temperature above 10°C for at least 2 hours and then completely inverted three times or placed in warm water at body temperature for 5 minutes and then completely inverted three times.

MAGNETIC RESONANCE (MR) CONDITIONAL REQUIREMENTS:

GCE medical combination valve itself is MR conditional. It means that this medical device can be safely used in an MR environment meeting the following conditions:

According to the test results the device can be considered safe for static magnetic fields up to 3 Tesla, with a maximum local static magnetic field of 40 mT (400 Gauss) or a spatial field gradient of $|\nabla B| < 0.5$ Tesla/meter.

Note: GCE cannot guarantee the MR conditional status of the whole gas package.

4. PERSONNEL INSTRUCTIONS

The medical device provider must ensure that all personnel handling the medical device are provided with the operating instructions & performance data. Gas source is not an integrated part of the medical device.

4.1. INTENDED USER

The intended user is any person, i. e. Healthcare professional or lay user*, who has read and understood the instructions for use and is able to operate the device correctly.

*The device is intended for lay users whose ability to operate the device is approved by the healthcare professional who prescribed the gas therapy.

If a training is considered necessary, it should be provided by the provider of the gas package (gas cylinder + integrated combination valve + associated equipment (cannula, breathing mask). The device shall be installed by a filling station and delivered to the end-user as ready to use the device.

 **Do not use the medical device without proper familiarization with the medical device and its safe operation as defined in this Instructions for use. Be aware of particular information and knowledge required for the gas in use.**

Note: In case of doubt regarding the use of the device as described by this Instructions for use, contact the medical device provider or manufacturer.

Note: In case of doubt regarding the gas therapy or in case of any undesired change in health condition, contact your healthcare professional.

5. MEDICAL DEVICE DESCRIPTION

A medical combination valve combines the function of shut-off valve of a high-pressure gas cylinder and pressure regulator for the use with medical gases. Gas from the cylinder is first controlled by the shut-off valve and then passed through the pressure regulator to the flow outlet or the pressure outlet, either of which can supply gas to a patient. The outlet pressure is fixed and downstream overpressure protected by a low-pressure relief valve.


EN

There are three basic alternatives:

- combination valve with a quick coupler outlet, the outlet pressure is fixed and the flow provided is on the demand of the downstream equipment, (there is no selectable flow outlet)
- combination valve with flow outlet, a selectable range of flows are controlled by the combination valve and can be changed by the flow control head, (there is no quick coupler outlet)
- the combination valve containing both quick coupler outlet and flow outlet.


Detailed information about the performance and technical data related to your device can be found in Appendix Nr. 1. Appendix Nr. 1 is in form of symbols. Explanations for the symbols are available in Chapter 12 - Glossary.

SHUT OFF VALVE STATUS INDICATOR:

Typical MediVital® valve configuration

MARKING DETAIL

Note: The content of the marking is always the same, the position on the valve may vary depending on the type.

A – INLET STEM

The combination valve is fitted to the gas cylinder by a threaded inlet stem. The inlet stem can be taper threaded or parallel threaded with different sizes depending on the cylinder size and material.

B – FILLING PORT

A filling port is provided for filling the gas cylinder at a filling station, it has no function to patient use. The filling port contains a non-return valve (NRV). The NRV prevents flow of gas without the use of specific filling adaptors which enables venting/evacuating gas from the cylinder and filling gas into the cylinder.

C – PRESSURE INDICATOR

The pressure indicator displays the amount of gas in the gas cylinder. The pressure indicator is of an active type which means it always indicates amount of gas in the gas cylinder whether the shut-off valve is opened or closed.

D – SHUT-OFF VALVE

The combination valve is provided with a shut-off valve to isolate the gas in the cylinder from the rest of the combination valve functions. It must be opened during cylinder filling and gas supply. Part of the handwheel is an indicator showing open/close status of the shut-off valve.

 The shut-off valve open/close status indicator is for guidance only. The shut-off valve may not be fully off when OFF status is showing. Fully closed status has to be checked by ensuring the shut-off valve is turned fully clockwise until a firm resistance to rotation is reached and gas supply to the outlets has stopped.

 The shut-off valve must not be used in the "partly ON" status because even though gas is supplied to the outlets, the flow can be limited due to an insufficiently opened shut-off valve.

E – RESIDUAL PRESSURE DEVICE

The combination valve is equipped with a residual pressure device (RPD) with a function to retain a minimum positive pressure in the gas cylinder to avoid contamination of the cylinder. During cylinder gas venting/evacuation through the filling port using a specific filling adapter the residual pressure device is bi-passed.

F, G – FLOW CONTROL HEAD “F” AND FLOW OUTLET “G” (OPTIONAL)

The combination valve may be fitted with a flow control head “F”. This function is used to supply gas flow rates (l/min) at atmospheric pressure directly to a patient through the flow outlet “G”, a range of selectable flow rates are provided. The flow outlet is connected with a cannula tube to provide gas to a patient either through a nasal cannula outlet or facemask. The flow rate selected is shown by the value identified by the triangular pointer. The flow outlet “G” is equipped with a hose fitting (hose nipple) or a threaded type (for related medical devices to be connected via threaded connection). Movement of the flow outlet ‘G’ is normal due to the method of fixing in the main body. It doesn’t indicate a faulty flow outlet.

H – PRESSURE OUTLET OR QUICK COUPLER (OPTIONAL)

The combination valve may be fitted with a pressure outlet. The pressure outlet is supplied with gas direct from the low-pressure part of the combination valve and it is fitted with a gas specific medical quick coupling connector also called “quick coupler”. The user can connect another piece of equipment to this outlet with a gas specific male probe. The quick coupler self-seals when the male probe is disconnected. This outlet is for supplying gas to other medical devices.

I - LOW PRESSURE RELIEF VALVE

The pressure relief valve provides protection from overpressure of the low-pressure part of the combination valve and connected medical devices. The low pressure relief valve vents gas if an overpressure scenario occurs until the gas pressure returns to normal conditions when it shall deactivate.

J - EXCESS FLOW DEVICE OR ANTI-DUST TUBE (OPTIONAL)

An excess flow device provides additional protection if the main part of the combination valve is sheared off above the inlet stem. The excess flow device activates in such situations to reduce the flow out of the inlet stem.

An anti-dust tube may be provided, its main function is to avoid contaminations within the gas cylinder entering the combination valve. It provides this function by collecting the gas at a more centralised position than the inlet stem provides and by a reduced aperture. The excess flow device also provides a contamination protection function. An anti-dust tube does not provide the same flow limiting function as an excess flow device.

PRESSURE RELIEF DEVICE OF HIGH PRESSURE PART - BURSTING DISC (OPTIONAL)

A high pressure relief device is intended to protect the cylinder and high-pressure part of the combination valve against damage caused by increased cylinder pressure. If the bursting disc has been activated, it does not reseal, and the combination valve must be taken out of service for repair (see Chapter 10).

GUARD (OPTIONAL)

A guard may be supplied as an integral part of the combination valve. The guard consists of two halves assembled together to the gas cylinder.

Note: Colour of the medical device (especially guard, flow control head and shut-off valve) does not have to match the gas colour coding.

 For manipulation use only the guard handle. Do not use other parts of the valve as handle.

6. INSTALLATION

See appendix Nr 4 – Valve assembly and filling instruction & appendix Nr 3 – Post Filling Checks.

7. OPERATION

EN

The medical gas is not an integrated part of the medical device. Gas therapy must be prescribed by a healthcare professional only.

7.1. OPERATIONS PERFORMED BY THE USER (BEFORE USE)

Contact a relevant filling station or manufacturer if you need help with the maintenance or unexpected operation of the combination valve.

7.1.1. VISUAL INSPECTION BEFORE USE:

- Check the combination valve for damage (incl. label, marking and printed marking). If it shows signs of external damages, remove the medical device from service and suitably identify its status.
- Check the combination valve for contamination. If needed apply the cleaning process described in this instructions for use.
- Check that the cylinder gas pressure indicator indicates sufficient pressure. If it indicates in the red zone, return the cylinder with combination valve back for filling.

7.1.2. LEAKAGE CHECK AND FUNCTIONAL TEST BEFORE USE:

- If applicable remove the protection elements (e.g. caps, plugs, covers). Keep them in a safe place for reuse during transport or storage.
- Set flow control head (if any) on the ZERO position - ensure the flow control head engages correctly.
- Slowly open the shut-off valve handwheel (anticlockwise) until fully open - after approx. 1 turn.
- Check for leakage by listening (leakage would be heard as the characteristic hiss of flowing gas).
- Check that there is a gas flow at each flow control head set position in both clockwise and anticlockwise turning direction (for instance by sound or checking presence of bubbles in a humidifier).
- Close the shut-off valve (clockwise). Do not use excessive torque (max. closing torque is 5 Nm).
- Reset the flow control head to ZERO position and ensure the flow control head engages correctly.
- For a combination valve fitted with a pressure outlet, ensure it is in working condition by connecting and disconnecting quick-coupler probe.

If any leakage is detected, use the procedure in chapter 7.3 and return the medical device to the provider.

Sudden opening could result in a danger of fire or explosion arising from oxygen pressure shocks.

7.2. USE OF THE COMBINATION VALVE

7.2.1. USE OF THE COMBINATION VALVE FLOW OUTLET AND SETTING OF FLOW:

- Ensure the flow control head is in position "O".
- Ensure that the hose/cannula is compatible with the flow outlet features and is well secured.
- Connect associated equipment to the flow outlet.
- Slowly open shut-off valve handwheel (anticlockwise) until fully open - after approx. 1 turn (see indication in picture chapter 5).
- Set flow head to required flow position which corresponds to value pointed by the triangular pointer. Ensure flow control head engages correctly.

Before connecting any associated equipment to the flow outlet make sure that the patient is not connected.

-  Sudden opening could result in a danger of fire or explosion arising from oxygen pressure shocks. Insufficient opening of the shut-off valve could reduce actual flow delivered.
-  Always ensure that the flow control head has engaged and not placed between two settings otherwise the combination valve will not deliver the correct flow of medical gas.
-  Common variants of the flow control head can have an “end stop” in between the maximum flow position and the zero position. Do not try to apply excessive force on the flow control head when it stops in the maximum flow position or in the zero position.
-  Do not use the flow outlet for driving any medical equipment.

7.2.2. USE OF COMBINATION VALVE PRESSURE OUTLET

- Ensure the flow control head is in position “0” (if included).
 - Ensure that the associated equipment IS NOT connected to the pressure outlet.
 - Slowly open the shut-off valve handwheel (anticlockwise) until fully open - after approx. 1 turn (see indication in picture chapter 5).
 - Ensure the medical quick coupler of associated equipment is compatible with the pressure outlet feature.
 - Connect the medical quick coupler of associated equipment into the pressure outlet.
 - Ensure that quick coupler of associated equipment is well secured.
-  Before connecting any associated equipment to the pressure outlet make sure that the patient is not connected and the associated equipment outlet is secured.
 -  If the pressure outlet is to be connected to a medical device that requires a high gas flow (for instance a pulmonary ventilator that requires a gas flow of 100 l/min at the minimum pressure 2.8 bar), compare the required flow of connecting the medical device with the pressure and flow characteristics of the combination valve stated in appendix Nr. 1. To assure sufficient performance (pressure and flow characteristics of the combination valve) the medical device should not be used if the pressure indicator enters the red zone.
 -  Sudden opening could result in a danger of fire or explosion arising from oxygen pressure shocks. Insufficient opening of the shut-off valve could reduce the actual flow delivered.
 -  If a pressure outlet and flow outlet are used simultaneously, GCE does not guarantee the outlet parameters as specified in Appendix Nr. 1 - Flow curve.

7.3. AFTER USE

- Close the shut-off valve (clockwise). Do not use excessive torque(max. closing torque is 5 Nm).
- Vent pressure from the connected devices.
- Disconnect all connected devices from user outlets.
- Set flow control head on "0" (if included).
- Refit the protection elements back (if present). Ensure that the protection elements are clean and undamaged.

8. ASSOCIATED EQUIPMENT

ASSOCIATED EQUIPMENT CONNECTABLE TO FLOW OUTLET:

- Hose connected with mask, cannula or humidifier.

ASSOCIATED EQUIPMENT CONNECTABLE TO PRESSURE OUTLET:

- Low pressure hose (working pressure >10 bar), flowmeters, Venturi suction ejectors, lung ventilators.

OTHER USER ASSOCIATED EQUIPMENT:

- bed hanger, humidifier holder.

ASSOCIATED EQUIPMENT FOR FILLING STATIONS:

- Filling adaptor.

-  Before connecting any associated equipment or medical device to the combination valve, always check that they are fully compatible with connection features & performances of the medical device.
-  During use ensure that the connected hose is not mechanically damaged, throttled, or broken, which would reduce or break the supply of gas to the patient.

9. CLEANING

The medical device is not sterile and does not require cleaning for contamination control. The purpose of cleaning is solely to remove contamination on the external surfaces (e.g. physical dirt) that could potentially cause safety issues with the medical device.

Remove dirt with a soft cloth dampened in oil free soap water & rinsed with clean water. Cleaning can be carried out with an alcohol-based solution (with dampened wipes). If other cleaning solutions are used, check that they are not abrasive and they are compatible with the medical device materials (including labels) and gas (convenient cleaning solution - i.e. Meliseptol).

-  If the cleaning instruction is not strictly followed, there is the fading of the printed marking risk which may mislead the user.
-  Do not use cleaning solutions containing ammonia!
-  Do not immerse in water or any liquid.
-  Do not expose to high temperature (such as autoclave).
-  To apply the cleaning solution do not spray it as the spray may enter into the inner parts of the combination valve and cause contamination or damage.
-  Do not use pressure washer as it could damage or contaminate the combination valve.
-  If the inner parts of the combination valve have been contaminated do not continue to use the combination valve under any circumstances. It must be withdrawn from service.

10. MAINTENANCE

10.1. REPAIR

Repairs and service can be only done by a GCE certified person who also holds all necessary certificates in compliance with national standards for mounting and repair of dedicated gas devices. For information about a service in your area please contact GCE or distributor of GCE medical device. Usually the pressure regulators integrated with cylinder valves can be repaired while fitted on the cylinder.

Repair that does not need to be done by certified personnel include exchange of the below mentioned components/spare parts:

- guard,
- labels,
- protective covers and removable hose adaptors.

-  All labels on the equipment must be kept in good, legible condition by the provider and the user during the entire product lifetime.
-  Use only original GCE components/spare parts!

Any medical device sent back to GCE (or GCE authorised centre) for repair or maintenance shall be properly packaged to prevent contamination or damage during storage, transportation and handling.

For medical device to be repaired, the fault short description or some reference to a claim nr. should be indicated.

10.2. SERVICE

There is no determined preventive check for combination valves during the entire life of the valve.

11. MEDICAL DEVICE LIFETIME

11.1. LIFETIME

Maximum lifetime (also known as service life) of the medical device is 15 years from the manufacturing date. At the end of the medical device's lifetime (15 years maximum), the medical device must be withdrawn from service.

11.2. SERIAL NUMBER AND DATE OF MEDICAL DEVICE

Form of nine-digit serial number stamped on the medical device is following:

YY MM XXXXX

YY: year of medical device

MM: month of medical device

XXXXX: sequence number

Example: serial number 220300001 shows the medical device produced in March 2022, with sequence number 1.

11.3. MARKING UDI CODING IMPLEMENTATION

There are 4 marked identifiers on the medical device surface:

14 numbers with prefix (01): GTIN - Global Trade Identification Number (identification of the selector type)

10 numbers with prefix (10): batch number (LOT)

9 numbers with prefix (21): serial number in format YY MM XXXXX

6 numbers with prefix (11): the medical device's production date in format YYMMDD

Example: marking number (11) 220308 shows the medical device with the production date of 8th March 2022.

As an integral part of UDI marking is a Data Matrix, which consists of the same data.

11.4. DISPOSAL OF WASTE MANAGEMENT

The owner of the device shall prevent the reuse of the medical device and handle the medical device in compliance with "Directive of European Parliament and Council 2008/98/EC on waste ". Contact your local GCE representative for further details, before returning or removing the medical device. All components should be disposed of according to local environmental laws applicable in the country of disposal. If required, a list of component materials is available from GCE.

11.5. REACH AND ROHS

In accordance with Article 33 of REACH GCE, s.r.o. as the legal manufacturer shall inform all customers if any materials exist containing 0.1% or more of substances included in the list of Substance of Very High Concern (SVHC). The most commonly used brass alloys for valve bodies and other brass components contain 2-3% of lead (Pb), EC no. 231-100-4, CAS no. 7439-92-1.

The lead will not be released to the gas or surrounding environment during normal use.

After the end of life the medical device shall be scrapped by an authorized metal recycler to ensure efficient material handling with minimal impact on the environment and health. To date we have no information that indicates that other materials containing SVHC of concentrations exceeding 0.1% are included in any GCE medical device.

11.6. REPORTING OF A SERIOUS INCIDENT

If any suspected serious incident has been occurred in relation to this medical device, notify the manufacturer, e-mail: adverse_events@gcegroup.com and the competent authority of the State in which the user and/or patient is established.

By reporting a suspected serious incident, you can help obtain more information about the safety of this medical device.

12. GLOSSARY

EN

	Consult instructions for use		Suitable for Home care use
	Caution		Suitable for Hospital care use
	Keep away from heat and flammable materials		Suitable for Emergency care use
	Manufacturer		Temperature limit
REF	Catalogue number	REF 1	Customer number
LOT	Batch code		Fragile, handle with care
	Inlet parameter		Outlet parameter
	Keep away from oil and grease!		Keep dry
MD	Medical Device		Weight of medical device
P₁	Inlet pressure		Contains hazardous substances
P₂	Outlet pressure	BD	Bursting disc
P₄	Maximal outlet pressure (closing pressure)		Date of manufacture
CH REP	Authorised representative for Switzerland	UDI	Unique Device Identification
		SN	Serial number

13. WARRANTY

The Standard Warranty period is 2 years from date of receipt by the GCE Customer (or if this is not known 2 years from time of the medical device manufacture shown on the medical device). The standard warranty is only valid for medical devices handled according to Instructions for use (IFU) and general industry good practice and standards.

APPENDIX:

Nr. 1 – Technical and performance data

Nr. 2 – Quick-coupler features and connecting / disconnecting procedure

Nr. 3 – Post Filling Checks

Nr. 4 – Valve assembly and filling instructions

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© GCE, s.r.o.